axiostat[®] stops bleeding, instantly

100% chitosan HAEMOSTATIC dressing

dedicated to Stop bleeding







Axiostat® is a sterile, single use, non-absorbable haemostatic dressing intended to be used for temporary control of moderate to severe bleeding wounds.

chitosan dressing made from pure grade chitosan extracted from shellfish. Chitosan is a linear polysaccharide consisting of glucosamine and N-acetyl glucosamine chains and is highly bio-active. Cationic nature of chitosan makes

Axiostat® is manufactured in **ISO 13485:2016** and **GMP** certified facility.

it mucoadhesive and thereby stops bleeding instantly.

Axiostat® is a 100%

Axiostat® is **India's first**

USFDA approved Haemostat

- Effective on patients on blood thinners
- Active on all sides
- Hypoallergenic
 no known side-effects
- O Can be cut, folded, stuffed into deep wounds
- O Gamma Sterilized
- Biocompatible
- O Easy to Apply & Remove

CLINICAL **DATA**

Haemostasis time in minutes

2 4 6 8 10 12 14

Standard
Gauze

13.08

Axiostat® 2

Clinical Study on Trauma patients Published in Indian Journal of Emergency Medicine IJEM VOL.2 No.2 JULY-DEC 2016

HOW IT WORKS



+ ly charged Axiostat®

Axiostat® reacts with blood cells and creates a robust seal



For detailed Cases/Clinical Studies, Visit: www.axiobio.com/case-studies or Scan this Code

AXIOSTAT® APPLICATIONS











L88 - 8x8cm | M85 - 8x5cm S55 - 5x5cm | V55 - 5x5cm



VASCULAR

V55 - 5x5cm V35 - 3.5x3.5cm

SURGICAL



All ectomy cases



General Surgery

Suture line anastomosis | Arterial and venous spurting Cardiothoracic & Vascular Surgeries (CTVS) | GI Surgery Post resection tumour bed bleeding | Capillary oozing Donor site bleeding | Orthopedic & Traumatic Surgeries



Interventional Radiology procedures



Interventional Cardiology procedures (Femoral & Radial)

Any type of Vascular interventional procedures



MIL88 - 8x8cm



Gunshot wounds



Puncture or Stab wounds

Blast injuries | Scalp and Limb cuts | Avulsions



N22 - 2x2cm

HAEMODIALYSIS



procedure



Haemodialysis

Vascular access procedures | Intravenous catheters Arterial/Central Line Removal



L88 - 8x8cm | M85 - 8x5cm S55 - 5x5cm | N14 - 1x4cm

N18 - 1x8cm





Head & Neck Trauma



Arterial bleeding

Venous Bleeding | Scalp Injuries | Cuts & Lacerations Limb injuries



D11 - 1x1cm | N22 - 2x2cm N14 - 1x4cm | V35 - 3.5x3.5cm



Maxillofacial trauma



Oral Surgeries

Extraction sockets | Oral biopsy | Flap surgery Gingival bleeding | Multiple tooth extractions Commando surgery & Tumour removal | Other oral procedures

___bleeding __use axiostat

ORDERING INFORMATION

PRODUCT CODE	PRODUCT SIZE	PACKAGING
L88	8cm x 8cm	Box of 10
M85	8cm x 5cm	Box of 10
S55	5cm x 5cm	Box of 10
N18	1cm x 8cm	Box of 20, 40
N14	1cm x 4cm	Box of 40, 80
MIL88	8cm x 8cm	Box of 10
V55	5cm x 5cm	Box of 10
V35	3.5cm x 3.5cm	Box of 15
N22	2cm x 2cm	Box of 20, 80
D11	1cm x 1cm	Box of 20, 80



AXIO BIOSOLUTIONS PRIVATE LIMITED

Plot 18, Gujarat Pharma Techno Park, Sari Matoda, Sanand Taluka, Ahmedabad - 382220, Gujarat, India Customer care: +91 8860786067, Email: info@axiobio.com | www.axiobio.com | twitter.com/axiobio facebook.com/axiobio